



# ELK Biotechnology

TAP-Tag Mouse mAb


Catalog NO.: EM1127

For research use only.

## Overview

---

Product name	TAP-Tag Mouse Monoclonal antibody
Source	Mouse
Applications	<b>WB IP</b>
Species reactivity	<b>N/A</b>
Recommended dilutions	<b>WesternBlot:1/5000-10000</b> <b>Immunoprecipitation:1/200</b> <b>NOTE: Optimal dilutions should be determined by the end user.</b>
Immunogen	Recombinant Protein
Species	N/A
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20° C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1 mg/ml
Observed band	<b>N/A</b>
GeneID (Human)	N/A
Human Swiss-Prot No.	N/A
Cellular localization	N/A
Alternative Names	tandem affinity purification
Background	Epitope mapping within the calmodulin binding peptide (CBP) region of the TAP (tandem affinity purification) tag of human.



Western blot analysis of TAP recombinant protein with EM1127 diluted at:5000.