

SRT-Tag Mouse mAb Catalog NO.: EM1154 For research use only.

Overview

Product name SRT-Tag Mouse Monoclonal antibody

Source Mouse

Applications WB

Species reactivity N/A

Recommended dilutions WesternBlot:1/5000

NOTE: Optimal dilutions should be determined by the end user.

Immunogen Recombinant Protein

Species N/A

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Store at -20° C. Avoid repeated freeze-thaw cycles.

Isotype IgG1

Clonality Monoclonal

Concentration 1 mg/ml

Observed band N/A

GenelD (Human) N/A


Human Swiss-Prot No. N/A

Cellular localization N/A

Alternative Names SRT Tag antibody

Background The SRT-tag is a peptide epitope comprised of 0 amino acids(TFIGAIATDT)

derived from the crystalline surface layer protein of Rickettsia typhi1.



Western blot analysis of SRT-Recombinant protein with EM1154 diluted at): $5000\ 2$): 10000