

S1-Tag Mouse mAb Catalog NO.: EM1195 For research use only.

Overview

Product name S1-Tag Mouse Monoclonal antibody

Source Mouse

Applications WB

Species reactivity ALL

Recommended dilutions WesternBlot:1/5000-10000

NOTE: Optimal dilutions should be determined by the end user.

Immunogen Synthetic Peptide

Species N/A

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Store at -20° C. Avoid repeated freeze-thaw cycles.

lsotype lgG1

Clonality Monoclonal

Concentration 1 mg/ml

Observed band N/A

GenelD (Human) N/A

Human Swiss-Prot No. N/A

Cellular localization N/A

Alternative Names N/A

Background S1 tag is an epitope tag composed of a nine-residue peptide NANNPDWDF

derived from the hepatitis B virus preS1 region.



Western blot analysis of Recombinant S1-Tag Protein with S1-Tag Mouse mAb EM1195 diluted at: 5000.