



## **ELK Biotechnology**

### **CY3-conjugated Goat Anti-Rabbit IgG**

**Catalog NO. : SA005**

**For research use only.**

#### **Overview**

---

<b>Product name</b>	CY3-conjugated Goat Anti-Rabbit IgG
<b>Recommended dilutions</b>	<b>1:20-1:100 for most applications</b> <b>NOTE: Optimal dilutions should be determined by the end user.</b>
<b>Size</b>	100µL
<b>Storage</b>	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C for one year. Avoid repeated freeze-thaw cycles.
<b>Note</b>	This product is for research use only, not for diagnostic or therapeutic use.

---

<b>Fluorophore</b>	Cy3
<b>Absorption Peak</b>	550nm
<b>Emission Peak</b>	570nm